

microBLOCK

A breakthrough in size and economy, our microBLOCK Mini Dry Bath is truly the first personal block incubator.

Whether your customers are incubating PCR strips or tube sizes up to 50ml, the compact microBLOCK Mini has the smallest footprint of any digital dry bath on the market.

Includes: Microblock dry bath and power adapter



KEY FEATURES

- Exchangeable blocks, for tubes 0.2 to 50ml
- Clear cover ensures temp. uniformity
- Digital temperature control
- Outstanding heating rate
- Aluminium blocks



A-MINI



B-MINI



C-MINI/D-MINI



E-MINI



F-MINI



G-MINI

TECHNICAL SPECIFICATIONS

Unit Dimensions	12 x 15.2 x 11.2cm
Weight	850 g
Timer	0 - 999 min
Input Voltage	24V
Operating Temperature	5°C - 35°C
Temperature Range	Ambient +5°C to 100°C
Uniformity	+/- 0.5°C
Temperature Increment	0.1°C
Block Dimension (W x L x H)	71 x 47 x 32 mm

stirring water baths

A powerful magnetic stirring mechanism combined with high wattage heating allows each stirring water bath to maintain temperatures to a maximum 99°C

Available in 10 and 20 litre bath capacities, these water baths comprise as many as 3 stirrers for a maximum stirring speed of 1500rpm. Each bath includes a highly visible front-panel LCD, reproducible microprocessor control of temperature within 0.1°C increments, a corrosion resistant stainless steel interior and automatic alarm and safety shutdown mechanism.



KEY FEATURES

- Powerful magnetic stirring mechanism
- Stirring speed up to 1500rpm
- Available in 10 and 20 litre
- Reproducible temperature control within 0.1°C

TECHNICAL SPECIFICATIONS

SWB-	10L-1	10L-2	20L-1	20L-3
Stirrers	1	2	1	3
Capacity (approx.)	10 L	10 L	20 L	20 L
Internal Dimensions	24x30x15cm		30x50x15cm	
Temperature	5°C above ambient to 99°C			
Heating Power	600 W	600 W	800 W	800 W
Stirring Speed	400 - 1500 rpm			
Timer	up to 99hr 59min, continuous			
Temperature	0.1°C			
Safety	warning indicator on screen, with alarm and automatic shut down			

ORDERING INFORMATION

MBDB-01* microBlock Digital Dry Bath with block lifter (Blocks sold separately)

A-MINI Block, for 0.2ml tubes (PCR Strip Tube), 40 wells

B-MINI Block, for 0.5ml tubes, 24 wells

C-MINI Block, for 15ml tubes, 15 wells

SWB-10L-1* Stirring Water Bath 10L with 1 built-in stirrer, includes lid

SWB-10L-2* Stirring Water Bath 10L with 2 built-in stirrers, includes lid

SWB-LID10 Transparent lid for 10L stirring water bath

* For 110V units, add \$ to order code

D-MINI Block, for 2ml tubes, 15 wells

E-MINI Block, for 12.5 ml cuvettes, 8 wells

F-MINI Block, for 15ml tubes, 4 wells

G-MINI Block, for 50ml tubes, 2 wells

SWB-20L-1* Stirring Water Bath 20L with 1 built-in stirrer, includes lid

SWB-20L-3* Stirring Water Bath 20L with 3 built-in stirrers, includes lid

SWB-LID20 Transparent lid for 20L stirring water bath